

## Comprehensive Epilepsy Panel

### Clinical Features

Epilepsy is defined by the occurrence of at least two unprovoked seizures occurring more than 24 hours apart. It is a common neurological disorder that affects at least 0.8% of the population. The International League against Epilepsy (ILAE) classifies seizures into two main categories.<sup>1</sup> Generalized epileptic seizures originate in and rapidly engage both cerebral hemispheres. Tonic-clonic, absence, myoclonic, clonic, tonic, and atonic seizures are common types of generalized seizures. Focal seizures originate from neuronal networks within a single hemisphere. Focal seizures are classified as focal aware (previously simple partial seizures), which do not result in an alteration of consciousness, and focal impaired (previously complex partial seizures), which cause a change in behavior or consciousness. Some types of seizures, such as infantile spasms, do not fit into either category and remain unclassified.

Seizures can be self-limiting or controlled by standard therapeutic treatments in some cases; however, individuals with epileptic encephalopathy have severe seizures that are refractory to treatment, leading to cognitive and behavioral impairment secondary to the epileptic activity. Epilepsy may be an isolated neurological symptom, or it may occur in association with other neurological symptoms or medical problems.<sup>2</sup> Some individuals with epilepsy are diagnosed with an electroclinical syndrome such as West or Ohtahara syndrome based on the presence of characteristic EEG findings and the clinical and family history.<sup>1</sup>

### Genetics

Epilepsy can be caused by genetic disorders, metabolic diseases, trauma, infection, and structural brain abnormalities, although the cause is not known in many cases. A genetic etiology underlies epilepsy in approximately 40% of individuals.<sup>3</sup> Genes have been identified that cause both generalized seizures and focal seizures, as well as unclassified epilepsy types such as infantile spasms. The genetic etiology of idiopathic generalized epilepsy (IGE) is frequently complex because it is due to a combination of multiple genetic factors that each confer a small risk for epilepsy and may be modified by environmental influences.<sup>3</sup> Currently, approximately 2% of patients with IGE harbor an identifiable variant in a single gene associated with Mendelian inheritance of epilepsy.<sup>4</sup> However, the percentage of patients with Mendelian epilepsy is higher for specific epilepsy types such as infantile spasms, benign familial neonatal and neonatal-infantile seizures (BFNS and BFNIS), and others.<sup>3,5,6,7</sup> The inheritance pattern can be autosomal dominant, autosomal recessive, or X-linked. Pathogenic variants in a single gene may be associated with different types of seizures (clinical heterogeneity), and conversely, pathogenic variants in different genes can cause the same epilepsy phenotype (genetic heterogeneity).

The Comprehensive Epilepsy panel at GeneDx includes sequencing and deletion/duplication analysis of genes causing Mendelian forms of epilepsy. Many of these genes encode subunits of ion channels involved in stabilizing or propagating neuronal activity, including components of the voltage-gated sodium, potassium, and calcium channels and the ligand-gated gamma-aminobutyric (GABA) channel.<sup>5,7,8,9</sup> The panel also includes non-ion channel genes associated with a variety of neurotransmitter disorders, storage

# TEST INFORMATION SHEET



and other neurometabolic disorders, as well as genes causing syndromic forms of epilepsy, many of which are involved in transcriptional activation or repression.<sup>5,6,8,10,11,12</sup>

## Test Methods

Using genomic DNA from the submitted specimen, the complete coding regions and splice site junctions of the genes on this panel are enriched using a proprietary targeted capture system developed by GeneDx for next-generation sequencing with CNV calling (NGS-CNV). The enriched targets are simultaneously sequenced with paired-end reads on an Illumina platform. Bi-directional sequence reads are assembled and aligned to reference sequences based on NCBI RefSeq transcripts and human genome build GRCh37/UCSC hg19. After gene specific filtering, data are analyzed to identify sequence variants and most deletions and duplications involving coding exons. Alternative sequencing or copy number detection methods are used to analyze or confirm regions with inadequate sequence or copy number data by NGS. Reportable variants include pathogenic variants, likely pathogenic variants and variants of uncertain significance. Likely benign and benign variants, if present, are not routinely reported but are available upon request.

The technical sensitivity of sequencing is estimated to be >99% at detecting single nucleotide events. It will not reliably detect deletions greater than 20 base pairs, insertions or rearrangements greater than 10 base pairs, or low-level mosaicism. The copy number assessment methods used with this test cannot reliably detect copy number variants of less than 500 base pairs or mosaicism and cannot identify balanced chromosome aberrations. Assessment of exon-level copy number events is dependent on the inherent sequence properties of the targeted regions, including shared homology and exon size. For the *CHRNA7* and *MAGI2* gene(s), deletion/duplication analysis but not sequencing is performed. For the *ARG1*, *C12orf57*, *FGF12*, *FOXG1*, *PHGDH*, *PIGG*, *PIGT*, *SLC35A2*, *SLC6A8*, and *SNAP25* gene(s), sequencing but not deletion/duplication analysis is performed. For the *ARX*, *ASNS*, *CHRNA7*, *FOLR1*, *GAMT*, *KCNT1*, *NR2F1*, *PCDH19*, *SCN1B*, *SHANK3*, and *TUBB2A* genes, only whole gene deletions or duplications may be detected, and deletions/duplications involving the 3' end of the *TSC2* gene (exons 36-42) may not be identified. If indicated based on the patient's specific clinical features, multiplex ligation-dependent probe amplification (MLPA) of the *FOXG1* or *SLC6A8* genes is available as a separate test (test code 906).

## Clinical Sensitivity

The clinical sensitivity of sequencing and deletion/duplication analysis of the genes included in the Comprehensive Epilepsy Panel depends in part on the patient's clinical phenotype. In a prior study, 31% of individuals with infantile spasms who were tested using an epilepsy gene panel were found to harbor definitive pathogenic variant(s) to explain the phenotype.<sup>113</sup> Overall, 19% of patients with epilepsy and 17-24% of patients with developmental and epileptic encephalopathies have an identifiable genetic etiology.<sup>114,175</sup> Specific information about the diagnostic yield for each gene in selected populations is summarized in the table below.

Appendix: Genes Included on Comprehensive Epilepsy Panel  
**Diagnostic Yield of Epilepsy Panel Genes in Selected Populations**

Epilepsy Type	Gene	Protein	Inh	Diagnostic Yield in Selected Population(s)
Benign familial neonatal seizures (BFNS)	KCNQ2	Potassium voltage-gated channel subfamily KQT member 2	AD	>50% BFNS <sup>6</sup>
	KCNQ3	Potassium voltage-gated channel subfamily KQT member 3	AD	~7% BFNS <sup>6</sup>
Benign familial neonatal-infantile seizures (BFNIS)	SCN2A	Sodium channel protein type 2 subunit alpha	AD	Unknown <sup>6</sup>
Benign familial infantile seizures (BFIS)	PRRT2	Proline-rich transmembrane protein 2	AD	83% familial and 30% sporadic BFIS <sup>14,15,16</sup> , 62–96% of familial PKD/IC; 36% sporadic PKD/IC <sup>14,17,18</sup>
Familial infantile myoclonic epilepsy (FIME)	TBC1D24	TBC1 domain family member 24	AR	Unknown <sup>19,20</sup>
Developmental and epileptic encephalopathy (includes West, Ohtahara, and Lennox Gastaut and other electroclinical syndromes)	ALDH7A1	Alpha-aminoacidic semialdehyde dehydrogenase (antiquitin)	AR	>90% pyridoxine-responsive epilepsy <sup>32</sup>
	ALG13	ALG13, UDP-N-acetylglucosaminyltransferase subunit	XL	Rare <sup>92</sup>
	ANKRD11	Ankyrin repeat domain-containing protein	AD	Rare <sup>154</sup>
	ARHGEF9	Cdc42 guanine nucleotide exchange factor (GEF) 9	XL	Rare <sup>94</sup>
	ARX	Aristaless related homeobox	XL	5% males with infantile spasms <sup>6</sup>
	ATP1A3	Sodium/potassium-transporting ATPase subunit alpha-3	AD	Rare in early life epilepsy <sup>131</sup> , 78% alternating hemiplegia of childhood <sup>132</sup>
	ATP6AP2	Renin receptor	XL	Rare <sup>45</sup>
	BRAT1	BRCA1 associated ATM activator 1	AR	Rare <sup>123</sup>
	CACNA1A	Voltage-dependent P/Q-type calcium channel subunit alpha-1A	AD	Rare <sup>95</sup>

# TEST INFORMATION SHEET



	<i>CACNA1E</i>	Voltage-dependent R-type calcium channel subunit alpha-1E	AD	Rare <sup>158</sup>
	<i>CACNA1G</i>	Voltage-dependent T-type calcium channel subunit alpha-1G	AD	Rare <sup>159</sup>
	<i>CASK</i>	Peripheral plasma membrane protein CASK	XL	Rare <sup>151</sup>
	<i>CDKL5</i>	Cyclin-dependent kinase-like 5	XL	10-17% infantile spasms <sup>6</sup>
	<i>CHD2</i>	Chromodomain helicase DNA binding protein 2	AD	1% of epileptic encephalopathy <sup>91</sup>
	<i>CLCN4</i>	Chloride voltage-gated channel 4	XL	Rare <sup>133</sup>
	<i>DNM1**</i>	Dynamin 1	AD	1% of epileptic encephalopathy <sup>92</sup>
	<i>DOCK7</i>	Dedicator of cytokinesis protein 7	AR	Rare <sup>164</sup>
	<i>EEF1A2</i>	Eukaryotic translation elongation factor 1 alpha 2	AD	Rare <sup>96, 97</sup>
	<i>FGF12</i>	Fibroblast growth factor 12	AD	Rare <sup>166</sup>
	<i>FRRS1L</i>	DOMON domain-containing protein FRRS1L	AR	Rare <sup>135</sup>
	<i>GABBR2</i>	Gamma-aminobutyric acid type B receptor subunit 2	AD	Rare <sup>118</sup>
	<i>GABRA1</i>	Gamma-aminobutyric acid receptor subunit alpha-1	AD	Rare <sup>88</sup>
	<i>GABRB2</i>	Gamma-aminobutyric acid (GABA) A receptor, beta 2	AD	Rare <sup>112</sup>
	<i>GABRB3</i>	Gamma-aminobutyric acid (GABA) A receptor, beta 3	AD	Rare <sup>95</sup>
	<i>GABRG2</i>	Gamma-aminobutyric acid receptor subunit gamma-2	AD	Rare <sup>6</sup>
	<i>GNAO1</i>	Guanine nucleotide binding protein (G protein), alpha activating activity polypeptide O	AD	Rare in epileptic encephalopathy and early-onset hyperkinetic movement disorders <sup>136</sup>
	<i>GRIN1</i>	Glutamate receptor, ionotropic, NMDA 1	AD	Rare <sup>94</sup>
	<i>GRIN2A</i>	Glutamate receptor ionotropic, NMDA 2A	AD	Rare <sup>83</sup>
	<i>GRIN2B</i>	Glutamate receptor ionotropic, NMDA 2B	AD	Rare <sup>90</sup>
	<i>HCN1</i>	Potassium/sodium hyperpolarization-activated	AD	Rare <sup>167</sup>

# TEST INFORMATION SHEET



	cyclic nucleotide-gated channel 1		
<i>IQSEC2</i>	IQ motif and sec7 domain 7	XL	Rare <sup>95</sup>
<i>KCNA2</i>	Potassium voltage-gated channel subfamily A member 2	AD	Rare <sup>138</sup>
<i>KCNB1</i>	Potassium channel, voltage gated shab related subfamily B, member 1	AD	Rare <sup>98</sup>
<i>KCNQ2</i>	Potassium voltage-gated channel subfamily KQT member 2	AD	10% neonatal epileptic encephalopathy <sup>31</sup>
<i>KCNT1</i>	Potassium channel, sodium activated subfamily T, member 1	AD	35% MMPSI <sup>99-101</sup> ; Rare in other epileptic encephalopathies <sup>99</sup>
<i>MAGI2 (dels only)</i>	Membrane-associated guanylate kinase, WW and PDZ domain-containing protein 2	AD	Rare <sup>43</sup>
<i>MEF2C</i>	Myocyte-specific enhancer factor 2C	AD	2% epileptic encephalopathy <sup>36</sup>
<i>NR2F1</i>	Nuclear receptor subfamily 2, group F, member 1	AD	Rare <sup>103</sup>
<i>PCDH19</i>	Protocadherin-19	XL	2-14% females with infantile/childhood epilepsy <sup>26,27,28,29,30</sup>
<i>PIGA</i>	Phosphatidylinositol glycan anchor biosynthesis, class A	XL	Rare <sup>109</sup>
<i>PIGG</i>	GPI ethanolamine phosphate transferase 2	AR	Rare <sup>169</sup>
<i>PIGN</i>	Phosphatidylinositol glycan anchor biosynthesis, class N	AR	Rare <sup>145</sup>
<i>PIGO</i>	Phosphatidylinositol glycan anchor biosynthesis, class O	AR	Rare <sup>110</sup>
<i>PIGT</i>	GPI transamidase component PIG-T	AR	Rare <sup>170</sup>
<i>PIGV</i>	Phosphatidylinositol glycan anchor biosynthesis, class V	AR	Rare <sup>111</sup>
<i>PLCB1</i>	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase beta-1	AR	Rare <sup>152</sup>
<i>PNPO</i>	Pyridoxine-5'-phosphate oxidase	AR	Rare <sup>40</sup>

# TEST INFORMATION SHEET



	<i>POLG</i>	DNA polymerase subunit gamma-1	AR	63–87% Alpers syndrome <sup>33,34,35</sup> ; 4–5% infantile/childhood epileptic encephalopathy <sup>34</sup>
	<i>QARS</i>	Glutaminyl-tRNA synthetase	AR	Rare <sup>104</sup>
	<i>SCN1A</i>	Sodium channel protein type 1 alpha	AD	70–80% Dravet syndrome <sup>6</sup> ; 20–24% early-onset cryptic epilepsy <sup>24,25</sup>
	<i>SCN2A</i>	Sodium channel protein type 2 alpha	AD	1–2% early-onset epileptic encephalopathy <sup>37,38</sup>
	<i>SCN8A</i>	Sodium channel protein type 8 subunit alpha	AD	Unknown <sup>39</sup>
	<i>SLC13A5</i>	Solute carrier family 13 (sodium-dependent citrate transporter), member 5	AR	Rare <sup>108</sup>
	<i>SLC25A22</i>	Mitochondrial glutamate carrier 1	AR	Rare <sup>44</sup>
	<i>SLC2A1</i>	Solute carrier family 2, facilitated glucose transporter member 1	AD	91% GLUT1 deficiency <sup>3</sup> ; ~10% early-onset absence epilepsy <sup>6</sup>
	<i>SLC35A2</i>	UDP-galactose translocator	XL	Rare <sup>172</sup>
	<i>SLC6A1</i>	Solute carrier family 6 member 1	AD	4% myoclonic-astatic epilepsy (MAE) <sup>149</sup>
	<i>SMC1A</i>	Structural maintenance of chromosomes protein 1A	AD	Rare <sup>151,152</sup>
	<i>SNAP25</i>	Synaptosomal-associated protein 25	AD	Rare <sup>173</sup>
	<i>SPATA5</i>	ATPase family protein 2 homolog	AR	Rare <sup>174</sup>
	<i>SPTAN1</i>	Alpha-II spectrin	AD	Rare <sup>41</sup>
	<i>STX1B</i>	Syntaxin 1B	AD	Rare in MAE and fever-associated epilepsy syndromes <sup>149</sup>
	<i>STXBPI</i>	Syntaxin binding protein 1	AD	35% Ohtahara syndrome <sup>6</sup> ; Unknown in Dravet syndrome <sup>88</sup>
	<i>SYNGAP1</i>	Ras/Rap GTPase-activating protein SynGAP	AD	Rare <sup>151</sup>
	<i>SZT2</i>	KICSTOR complex protein SZT2	AR	Rare <sup>149</sup>
	<i>TBL1XR1</i>	F-box-like/WD repeat-containing protein TBL1XR1	AD	Rare <sup>150</sup>
	<i>TSC1</i>	Hamartin	AD	2–4% infantile spasms <sup>21,22,23</sup>
	<i>TSC2</i>	Tuberin	AD	10–16% infantile spasms <sup>21,22,23</sup>
	<i>WWOX</i>	WW domain-containing oxidoreductase	AR	Rare

# TEST INFORMATION SHEET

GeneDx

Autosomal Dominant Focal Epilepsies	<i>CHRNA4</i>	Neuronal acetylcholine receptor alpha-4	AD	~10% autosomal dominant nocturnal frontal lobe epilepsy (ADNFLE) <sup>6</sup>
	<i>CHRN B2</i>	Neuronal acetylcholine receptor beta-2	AD	< 5% ADNFLE <sup>6</sup>
	<i>CHRNA2</i>	Neuronal acetylcholine receptor alpha-2	AD	Rare in ADNFLE <sup>6</sup>
	<i>DEPDC5</i>	DEP domain containing 5	AD	5–37% autosomal dominant focal epilepsies, with or without cortical dysplasia <sup>126</sup>
	<i>KCNT1</i>	Potassium channel, sodium-activated subfamily T, member 1	AD	< 5% ADNFLE <sup>99,102</sup>
	<i>LGII</i>	Leucine-rich glioma-inactivated protein 1	AD	~50% autosomal dominant partial epilepsy with auditory features (ADTLE) <sup>6</sup>
	<i>NPRL3</i>	GATOR complex protein NPRL3	AD	~2% autosomal dominant focal epilepsies, with or without cortical dysplasia <sup>127</sup>
Generalized Epilepsy with Febrile Seizures Plus (GEFS+)	<i>SCN1A</i>	Sodium channel protein type 1 alpha	AD	5–10% GEFS+ <sup>6</sup>
	<i>SCN1B</i>	Sodium channel subunit beta-1	AD	<5% GEFS+ <sup>6</sup>
	<i>GABRG2</i>	Gamma-aminobutyric acid receptor subunit gamma-2	AD	<1% GEFS+ <sup>6</sup>
	<i>SCN2A</i>	Sodium channel protein type 2 alpha	AD	Rare <sup>5</sup>
Epilepsy with paroxysmal dyskinesia	<i>PRRT2</i>	Proline-rich transmembrane protein 2	AD	62–96% of familial PKD/infantile convulsions (PKD/IC); 36% sporadic PKD/IC <sup>14,17,18</sup>
	<i>KCNMA1</i>	Potassium calcium-activated channel subfamily M alpha 1	AD/ AR	Rare <sup>140</sup>
Progressive Myoclonic Epilepsy	<i>EPM2A</i>	Laforin	AR	53% Lafora disease <sup>48</sup>
	<i>NHLRC1</i> ( <i>EPM2B</i> )	NHL repeat-containing protein 1 (malin)	AR	40% Lafora disease <sup>48</sup>
	<i>CSTB*</i> ( <i>EPM1</i> )	Cystatin-B	AR	Rare sequencing variants in Unverricht-Lundborg disease <sup>49</sup>
	<i>KCTD7</i> ( <i>EPM3</i> )	BTB/POZ domain-containing protein KCTD7	AR	Rare <sup>50</sup>
	<i>SCARB2</i> ( <i>EPM4</i> )	Lysosome membrane protein 2	AR	~7% progressive myoclonic epilepsy <sup>51</sup> ; Unknown in action

# TEST INFORMATION SHEET



				myoclonus-renal failure syndrome <sup>52</sup>
	<i>GOSR2</i>	Golgi SNAP receptor complex member 2	AR	Rare <sup>53,54</sup>
	<i>KCNC1</i>	Potassium voltage-gated channel subfamily C member 1	AD	Rare KCNC1-associated myoclonus epilepsy and ataxia (MEAK) <sup>115</sup>
	<i>FOLR1</i>	Folate receptor alpha	AR	Rare <sup>55</sup>
Rett/atypical Rett syndromes	<i>MECP2</i>	Methyl CpG binding protein 2	XL	88% females with Rett syndrome <sup>57</sup>
	<i>CDKL5</i>	Cyclin-dependent kinase-like 5	XL	2-8% females with atypical Rett syndrome <sup>58,59</sup>
	<i>CTNNB1</i>	Catenin beta 1	AD	Rare in Rett-like syndromes <sup>116</sup>
	<i>DDX3X</i>	DEAD-box helicase 3, X-linked	XL	Rare in Rett-like syndromes <sup>117</sup>
	<i>FOXG1**</i>	Forkhead box protein G1	AD	~1% Rett syndrome overall <sup>60</sup> ; 25% congenital variant of Rett <sup>61</sup>
	<i>GABBR2</i>	Gamma-aminobutyric acid type B receptor subunit 2	AD	Rare in epileptic encephalopathy and intellectual disability with Rett-like features <sup>118</sup>
	<i>KCNA2</i>	Potassium voltage-gated channel subfamily A member 2	AD	Rare in Rett-like syndromes <sup>119</sup>
	<i>MBD5</i>	Methyl-CpG-binding domain protein 5	AD	Rare in Rett-like syndromes <sup>62</sup>
	<i>MEF2C</i>	Myocyte-specific enhancer factor 2C	AD	Rare in Rett-like syndromes <sup>63</sup>
	<i>SATB2</i>	SATB homeobox 2	AD	Rare in Rett-like syndromes <sup>120</sup>
	<i>STXBP1</i>	Syntaxin binding protein 1	AD	Rare in Rett-like syndromes <sup>121</sup>
	<i>TBLIXR1</i>	Transducin (beta)-like 1 X-linked receptor 1	AD	Rare in Rett-like syndromes <sup>122</sup>
Angelman/ Angelman-like/ Pitt-Hopkins syndromes	<i>WDR45</i>	WD repeat domain 45	XL	Rare in Rett-like syndromes <sup>106,107</sup>
	<i>UBE3A</i>	Ubiquitin protein ligase E3A	AD	68% maternally inherited 15q11.2 deletion <sup>64</sup> ; 11% UBE3A sequencing variant Angelman syndrome <sup>64</sup>
	<i>SLC9A6</i>	Sodium/hydrogen exchanger 6	XL	6% Angelman-like syndrome <sup>65</sup>
	<i>MBD5</i>	Methyl-CpG-binding domain protein 5	AD	Rare in Angelman-like syndrome <sup>66</sup>
	<i>TCF4</i>	Transcription factor 4	AD	36% Pitt-Hopkins syndrome (PHS) <sup>67</sup> ; 2% Angelman syndrome <sup>67</sup>

# TEST INFORMATION SHEET



	<i>ATRX</i>	Transcriptional regulator ATRX	XL	95% males with alpha-thalassemia X-linked intellectual disability syndrome <sup>153</sup>
	<i>NRXN1</i>	Neurexin-1	AR	Rare in PHS <sup>68</sup>
	<i>CNTNAP2</i>	Contactin-associated protein-like 2	AR	Rare in PHS <sup>68</sup>
Creatine deficiency syndromes	<i>SLC6A8**</i>	Solute carrier family 6 (neurotransmitter transporter), member 8	XL	2% males with epilepsy and intellectual disability <sup>46</sup> ; 65% males with biochemical creatine deficiency <sup>47</sup>
	<i>GAMT**</i>	Guanidinoacetate N-methyltransferase	AR	Rare <sup>71</sup>
	<i>GATM</i>	Glycine amidinotransferase, mitochondrial	AR	Rare <sup>71</sup>
Neuronal Ceroid Lipofuscinoses (NCL)	<i>PPT1 (CLN1)</i>	Palmitoyl-protein thioesterase 1	AR	98% PPT1 deficiency <sup>72</sup>
	<i>TPP1 (CLN2)</i>	Tripeptidyl-peptidase 1	AR	95% TPP1 deficiency <sup>73</sup>
	<i>CLN3</i>	Battenin	AR	92% Juvenile NCL <sup>74</sup>
	<i>DNAJC5 (CLN4B)</i>	DnaJ homolog subfamily C member 5	AD	25% Kufs disease <sup>75,76</sup>
	<i>CLN5</i>	Ceroid-lipofuscinoses neuronal protein 5	AR	94% Finnish late-infantile NCL <sup>77</sup> ; Otherwise rare
	<i>CLN6</i>	Ceroid-lipofuscinoses neuronal protein 6	AR	Rare <sup>77</sup>
	<i>MFSD8 (CLN7)</i>	Major facilitator superfamily domain-containing protein 8	AR	Rare <sup>77</sup>
	<i>CLN8</i>	Ceroid-lipofuscinoses neuronal protein 8	AR	100% Finnish Northern epilepsy <sup>77</sup> ; Otherwise rare
	<i>CTSD (CLN10)</i>	Cathepsin D	AR	Rare <sup>77</sup>
	<i>CTSF (CLN13)</i>	Cathepsin F	AR	Rare <sup>77</sup>
Metabolic disorders	<i>KCTD7 (CLN14)</i>	BTB/POZ domain-containing protein KCTD7	AR	Rare <sup>78</sup>
	<i>ADSL</i>	Adenylosuccinate lyase	AR	Rare (adenosuccinate lyase deficiency) <sup>79</sup>
	<i>ALDH5A1</i>	Aldehyde dehydrogenase 5 family member A1	AR	96% (SSADH deficiency) <sup>129</sup>
	<i>ARG1</i>	Arginase 1	AR	Rare (argininemia) <sup>155</sup>
	<i>ASNS</i>	Asparagine synthetase (glutamine-hydrolyzing)	AR	Rare (asparagine synthetase deficiency) <sup>130</sup>
	<i>FOLR1</i>	Folate receptor alpha	AR	Rare (cerebral folate deficiency) <sup>56</sup>

# TEST INFORMATION SHEET

GeneDx

	<i>GLDC</i>	Glycine decarboxylase	AR	70–75% glycine encephalopathy (nonketotic hyperglycinemia) <sup>124</sup>
	<i>NGLY1</i>	N-glycanase 1	AR	Rare (congenital disorder of deglycosylation) <sup>143</sup>
	<i>PHGDH</i>	D-3-phosphoglycerate dehydrogenase	AR	Rare (phosphoglycerate dehydrogenase deficiency / Neu-Laxova syndrome) <sup>168</sup>
	<i>SLC19A3</i>	Solute carrier family 19 member 3	AR	Rare (thiamine metabolism dysfunction disorders) <sup>125</sup>
Brain malformations	<i>ARX</i>	Aristaless related homeobox	XL	70–95% of X-linked lissencephaly with ambiguous genitalia (XLAG) <sup>156</sup>
	<i>CASK</i>	Peripheral plasma membrane protein CASK	XL	~96% females with microcephaly with pontine and cerebellar hypoplasia (MICPCH) <sup>151</sup>
	<i>CUL4B</i>	Cullin 4B	XL	Rare in cerebral malformations <sup>160</sup>
	<i>DCX</i>	Doublecortin	XL	Up to 100% XL lissencephaly; 10% of classic lissencephaly; 85% females and ~30% males with subcortical band heterotopia (SBH) <sup>161–163</sup>
	<i>FLNA</i>	Filamin A	XL	93% classic bilateral periventricular nodular heterotopia <sup>134</sup>
	<i>PAFAH1B1 (LIS1)</i>	Platelet-activating factor acetylhydrolase IB subunit alpha	AD	~40–65% of classic lissencephaly <sup>162</sup>
	<i>TUBB2A</i>	Tubulin, Beta-2A	AD	Rare <sup>70</sup>
EAST/SeSAME syndrome	<i>KCNJ10</i>	ATP-sensitive inward rectifier potassium channel 10	AR	Rare <sup>81,82</sup>
Epilepsy and hemiplegic migraine	<i>ATP1A2</i>	Sodium/potassium-transporting ATPase subunit alpha-2	AD	Unknown <sup>6</sup>
Idiopathic generalized epilepsy	<i>CHRNA7 (dels only)</i>	Neuronal acetylcholine receptor subunit alpha-7	AD	~1% patients with IGE <sup>56</sup>
Epilepsy associated with other neuro-developmental disorders	<i>ANKRD11</i>	Ankyrin repeat domain-containing protein	AD	83% patients with KBG syndrome <sup>154</sup>
	<i>C12orf57</i>	Protein C10	AR	Rare (Temptamy syndrome or TEMTYS) <sup>157</sup>
	<i>DYRK1A</i>	Dual-specificity tyrosine-(Y)-phosphorylation regulated kinase 1A	AD	Rare <sup>93</sup>

# TEST INFORMATION SHEET



	<i>EEF1A2</i>	Eukaryotic translation elongation factor 1 alpha 2	AD	Rare <sup>96</sup>
	<i>EHMT1</i>	Histone-lysine N-methyltransferase EHMT1	AD	~75% of patients with Kleefstra syndrome, including microdeletion of 9q34.3 <sup>165</sup>
	<i>GRIN2A</i>	Glutamate receptor ionotropic, NMDA 2A	AD	9-20% atypical Rolandic epilepsy and epilepsy-aphasia <sup>89</sup> ; ~2% patients with intellectual disability and epilepsy <sup>83,84</sup>
	<i>GRIN2B</i>	Glutamate receptor ionotropic, NMDA 2B	AD	Unknown <sup>83</sup>
	<i>KANSL1</i>	KAT8 regulatory NSL complex subunit 1	AD	Rare (Koolen-de Vries syndrome, including microdeletion of 17q21.31) <sup>85,86</sup>
	<i>HNRNPU</i>	HNRNPU antisense RNA 1	AD	Rare <sup>137</sup>
	<i>KCNH1</i>	Potassium voltage-gated channel subfamily H member 1	AD	Rare (Temple-Baraitser and Zimmermann-Laband syndromes) <sup>139</sup>
	<i>KDM6A</i>	Lysine demethylase 6A	XL	Rare (Kabuki syndrome) <sup>128</sup>
	<i>NEXMIF (KIAA2022)</i>	Neurite extension and migration factor	XL	Rare <sup>141</sup>
	<i>NALCN</i>	Sodium leak channel, non-selective	AD/ AR	Rare (CLIFAHDD and infantile hypotonia with psychomotor retardation and characteristic facies (IHPRF) <sup>142</sup>
	<i>PACSI</i>	Phosphofuran acid cluster sorting protein 1	AD	Rare <sup>144</sup>
	<i>PNKP</i>	Bifunctional polynucleotide phosphatase/kinase	AR	Rare <sup>86</sup>
	<i>PPP2R5D</i>	protein phosphatase 2 regulatory subunit B', delta	AD	Rare <sup>146</sup>
	<i>PURA</i>	Transcriptional activator protein Pur-alpha	AD	Rare (Including microdeletion of 5q31.3) <sup>147</sup>
	<i>SHANK3</i>	SH3 and multiple ankyrin repeat domains protein 3	AD	Rare (Phelan-McDermid syndrome, including microdeletion of 22q13.3) <sup>171</sup>
	<i>SMARCA2</i>	Probable global transcription activator SNF2L2	AD	Rare (Nicolaides-Baraitser and Coffin Siris syndromes) <sup>176</sup>

## References

1. Berg et al. (2010) Epilepsia 51: 676-685.
2. Pellock, JM (2004) Neurol (2004) 62:S17-S23.
3. Pong et al., (2011) Pediatr Neurol 44:317-327.
4. Weber et al., (2008) Dev Med Child Neurol 50:648-654.
5. Nicita et al., (2011) Seizure: Eur J Epilepsy doi:10.1016/j.seizure.2011.08.007
6. Ottman et al., (2010) Epilepsia 51:655-670.
7. Pal et al., (2010) Nat Rev Neurol 6:445-453.
8. Deprez et al., (2009) Neurol 72:273-281.
9. Ottman et al., (2010) Epilepsia 51:655-670.



# TEST INFORMATION SHEET



al. (2015) Ann. Neurol. 77 (4):675–83 (PMID: 25623524) 127. Weckhuysen et al. (2016) Epilepsia 57 (6):994–1003 (PMID: 27173016) 128. Bögershausen et al. (2016) Hum. Mutat. 37 (9):847–64 (PMID: 27302555) 129. Akaboshi et al. (2003) Human Mutation 22 (6):442–50 (PMID: 14635103) 130. Gupta et al. (2017) Metab Brain Dis 32 (6):1889–1900 (PMID: 28776279) 131. Paciorkowski et al. (2015) Epilepsia 56 (3):422–30 (PMID: 25656163) 132. Heinzen et al. (2012) Nature Genetics 44 (9):1030–4 (PMID: 22842232) 133. Palmer et al. (2018) Mol. Psychiatry 23 (2):222–230 (PMID: 27550844) 134. Chen MH, Walsh CA. FLNA-Related Periventricular Nodular Heterotopia. 2002 Oct 8 [Updated 2015 Sep 17]. In: Adam MP, Ardinger HH, Pagon RA, et al., editors. GeneReviews® [Internet]. Seattle (WA): University of Washington, Seattle; 1993–2018. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK1213/> 135. Madeo et al. (2016) Am. J. Hum. Genet. 98 (6):1249–55 (PMID: 27236917) 136. Danti et al. (2017) Neurol Genet 3 (2):e143 (PMID: 28357411) 137. Bramswig et al. (2017) Hum. Genet. 136 (7):821–834 (PMID: 28393272) 138. Masnada et al. (2017) Brain 140 (9):2337–2354 (PMID: 29050392) 139. Mégarbané et al. (2016) BMC Med. Genet. 17 (1):42 (PMID: 27282200) 140. Webster et al. (2017) Clin. Genet. 91 (5):756–763 (PMID: 27568816) 141. Webster et al. (2017) Clin. Genet. 91 (5):756–763 (PMID: 27568816) 142. Bend et al. (2016) Neurology 87 (11):1131–9 (PMID: 27558372) 143. Enns et al. (2014) Genet. Med. 16 (10):751–8 (PMID: 24651605) 144. Schuurs-Hoeijmakers et al. (2016) Am. J. Med. Genet. A 170 (3):670–5 (PMID: 26842493) 145. Fleming et al. (2016) Am. J. Med. Genet. A 170 (1):77–86 (PMID: 26394714) 146. Shang et al. (2016) Neurogenetics 17 (1):43–9 (PMID: 26576547) 147. Reijnders et al. (2018) J. Med. Genet. 55 (2):104–113 (PMID: 29097605) 148. Carvill et al. (2015) Am. J. Hum. Genet. 96 (5):808–15 (PMID: 25865495) 149. Tsuchida et al. (2018) Clin. Genet. 93 (2):266–274 (PMID: 28556953) 150. Saitsu et al. (2014) Journal Of Human Genetics 59 (10):581–3 (PMID: 25102098) 151. Moog U, Uyanik G, Kutsche K. CASK-Related Disorders. 2013 Nov 26. In: Adam MP, Ardinger HH, Pagon RA, et al., editors. GeneReviews® [Internet]. Seattle (WA): University of Washington, Seattle; 1993–2018. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK169825/> 151. Mignot et al. (2016) J. Med. Genet. 53 (8):511–22 (PMID: 26989088) 152. Ngoh et al. (2014) Dev Med Child Neurol 56 (11):1124–8 (PMID: 24684524) 153. Stevenson RE. Alpha-Thalassemia X-Linked Intellectual Disability Syndrome. 2000 Jun 19 [Updated 2014 Nov 6]. In: Adam MP, Ardinger HH, Pagon RA, et al., editors. GeneReviews® [Internet]. Seattle (WA): University of Washington, Seattle; 1993–2018. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK1449/>. 154. Low et al. (2016) Am. J. Med. Genet. A 170 (11):2835–2846 (PMID: 27667800) 155. Jain-Ghai et al. (2011) Molecular Genetics And Metabolism 104 (1–2):107–11 (PMID: 21802329) 156. Dobyns (2010) Epilepsia 51 Suppl 1:5–9 (PMID: 20331703). 157. Zahrani et al. (2013) American Journal Of Human Genetics 92 (3):387–91 (PMID: 23453665) 158. Helbig et al. (2018) Am. J. Hum. Genet. 103 (5):666–678 (PMID: 30343943) 159. Kunii et al. (2020) J Neurol Sci 416 :117047 (PMID: 32736238) 160. Vulto-van et al. (2015) Hum. Mutat. 36 (1):106–17 (PMID: 25385192) 161. Dobyns (2010) Epilepsia 51 Suppl 1:5–9 (PMID: 20331703) 162. Fry et al. (2014) American Journal Of Medical Genetics. Part C, Seminars In Medical Genetics 166C (2):198–210 (PMID: 24862549). 163. Bahi-Buisson et al. (2013) Brain : A Journal Of Neurology 136 (Pt 1):223–44 (PMID: 23365099). 164. Haberlandt et al. (2021) Mol Genet Genomic Med :e1607 (PMID: 33471954) 165. Willemsen et al. (2012) Mol Syndromol 2 (3–5):202–212 (PMID: 22670141) 166. Shi et al. (2017) Neurol Genet 3 (1):e133 (PMID: 28144627). 167. Marini et al. (2018) Brain 141 (11):3160–3178 (PMID: 30351409). 168. Acuna-Hidalgo et al. (2014) Am. J. Hum. Genet. 95(3):285–93 (PMID: 25152457). 169. Makrythanasis et al. (2016) Am. J. Hum. Genet. 98 (4):615–26 (PMID: 26996948) 170. Bayat et al. (2019) Genet. Med. 21 (10):2216–2223 (PMID: 30976099). 171. Li et al. (2018) Am. J. Med. Genet. A 176 (12):2668–2676 (PMID: 30537371). 172. Yates et al. (2018) Eur. J. Paediatr. Neurol. 22 (6):1095–1102 (PMID: 30194038). 173. Klöckner et al. (2020) Genet Med : (PMID: 33299146) 174. Kurata et al. (2016) Clin. Genet. 90 (5):437–444 (PMID: 27246907) 175. Sheidley et al. (2022) Epilepsia 63(2):375–387. 176. Sousa et al. (2014) Am J Med Genet C Semin Med Genet 166C (3):302–14 (PMID: 25169058)