

## Rett/Angelman Syndromes and Related Disorders Panel

### PANEL GENE LIST

ATRX, CDKL5, CNTNAP2, CTNNB1, DDX3X, DYRK1A, EHMT1, FOXG1\*, GABBR2, IQSEC2, KCNA2, MBD5, MECP2, MEF2C, NRXN1, PCDH19\*\*, SATB2, SHANK3\*\*, SLC9A6, STXBP1, TBL1XR1, TCF4, UBE3A\*\*\*, WDR45, ZEB2

\* This panel does not include deletion/duplication testing of FOXG1

\*\* Only whole gene deletion/duplications may be detected for the PCDH19 and SHANK3 genes

\*\*\* This panel includes MS-MLPA to evaluate for abnormal methylation of the UBE3A gene

### CLINICAL FEATURES

The Rett/Angelman Syndrome and Related Disorders Panel at GeneDx includes genes that cause disorders with overlapping clinical phenotypes including epilepsy, developmental delay and/or regression, movement disorders, and intellectual disability. A brief summary of some of the clinical disorders is included below.

Rett syndrome is a progressive neurodevelopmental disorder that primarily affects females. Classic Rett syndrome is characterized by apparently normal development in the first 6-18 months followed by an arrest in development and subsequent regression in language and motor skills. Frequent symptoms include loss of speech and purposeful hand use, stereotypic hand movements, ataxia, microcephaly, and seizures. Multiple forms of atypical (variant) Rett syndrome have been described, including early-onset Rett syndrome with seizures beginning before six months, forme fruste Rett syndrome with a milder and more slowly progressive course, and the preserved speech variant with a milder clinical phenotype including retained verbal abilities and hand use.<sup>1</sup> Congenital Rett syndrome is a severe congenital encephalopathy that affects males and females and is characterized by hypotonia, irritability, and unresponsiveness in the neonatal period with seizures, progressive microcephaly and developmental delay/regression noted in the first few months of life.<sup>2,3</sup>

Angelman syndrome (AS) is characterized by developmental delay with absent or significantly impaired speech, intellectual disability, ataxia, microcephaly, a characteristic EEG pattern, and a typical behavioral pattern including a happy personality with outbursts of laughter, hyperactivity, excitability, sleep problems, and hand-flapping movements.<sup>4,5</sup> Facial features often observed in individuals with AS include a prominent chin, a wide mouth with a protruding tongue. Angelman-like (Christianson) syndrome is characterized by intellectual disability, ataxia, severe speech and language impairment, epilepsy, and microcephaly in males.<sup>6,7,8</sup>

Pitt-Hopkins syndrome (PHS) is characterized by severe intellectual disability with absent or severely delayed speech that is often associated with breathing abnormalities, stereotypic movements, seizures, microcephaly and a characteristic facial gestalt consisting of coarse facial features, deep-set eyes, a broad or beaked nasal bridge, a large mouth with a tented upper lip and widely spaced teeth, and cup-shaped ears with a thick helix.<sup>9,10,11</sup> Ataxia, short stature, constipation, gastroesophageal reflux, strabismus, and nonspecific abnormalities on brain MRI have also been described.<sup>9,10,11</sup>

Mowat-Wilson syndrome (MWS) is characterized by severe expressive language delay, seizures, moderate to severe intellectual disability, and a typical facial gestalt consisting of uplifted earlobes with a central depression, broad and medially thick eyebrows, hypertelorism, downslanting palpebral fissures, an open-mouthed expression, and an elongated face with prognathism.<sup>12,13</sup> Hirschsprung disease (HSCR), congenital heart defects, microcephaly, ataxia, urogenital anomalies, short stature, structural eye anomalies, and hypoplasia/aplasia of the corpus callosum have also been described in more than half of patients with MWS.<sup>12</sup>

## GENETICS

The Rett/Angelman Syndrome and Related Disorders Panel includes sequencing and deletion/duplication analysis for the genes listed above, as well as methylation-specific multiple ligation probe amplification (MS-MLPA) analysis of the UBE3A gene to also evaluate for paternal UPD (3%-7% of AS) and imprinting errors (~3% of AS) causing Angelman syndrome. The disorders on this panel are inherited in an autosomal dominant, autosomal recessive, or X-linked manner. The genes on this panel are highly expressed in the developing nervous system and encode proteins that are involved in transcriptional activation or repression, are implicated in cell recognition and adhesion, or are involved in regulation of endosomal membranes and long-term potentiation.

## TEST METHODS

Using genomic DNA from the submitted specimen, the complete coding regions and splice site junctions of the genes on this panel are enriched using a proprietary targeted capture system developed by GeneDx for next-generation sequencing with CNV calling (NGS-CNV). The enriched targets are simultaneously sequenced with paired-end reads on an Illumina platform. Bi-directional sequence reads are assembled and aligned to reference sequences based on NCBI RefSeq transcripts and human genome build GRCh37/UCSC hg19. After gene specific filtering, data are analyzed to identify sequence variants and most deletions and duplications involving coding exons; however, technical limitations and inherent sequence properties effectively reduce this resolution for some genes. Alternative sequencing or copy number detection methods are used to analyze or confirm regions with inadequate sequence or copy number data by NGS. Additionally, MS-MLPA is performed to evaluate for abnormal methylation of the UBE3A gene. Reportable variants include pathogenic variants, likely pathogenic variants and variants of uncertain significance. Likely benign and benign variants, if present, are not routinely reported but are available upon request.

The technical sensitivity of sequencing is estimated to be > 99% at detecting single nucleotide events. It will not reliably detect deletions greater than 20 base pairs, insertions or rearrangements greater than 10 base pairs, or low-level mosaicism. The copy number assessment methods used with this test cannot reliably detect copy number variants of less than 500 base pairs or mosaicism and cannot identify balanced chromosome aberrations. Assessment of exon-level copy number events is dependent on the inherent sequence properties of the targeted regions, including shared homology and exon size. For the FOXP1 gene(s), sequencing but not deletion/duplication analysis, is performed. For the PCDH19 and SHANK3 genes, only whole gene deletions or duplications may be detected. If indicated, multiplex ligation-dependent probe amplification (MLPA) of the FOXP1 gene is available as a separate test (test code 906).

## CLINICAL SENSITIVITY

The clinical sensitivity of sequencing and deletion/duplication analysis of the genes included in the Rett/Angelman Syndrome and Related Disorders depends in part on the patient's clinical phenotype. Specific information about the diagnostic yield for each gene in selected populations is summarized in the table below.

Diagnostic Yield of Rett/Angelman Panel Genes in Selected Populations

Disorder	Gene	Protein	Inh	Diagnostic Yield in Selected Population(s)
Rett / atypical Rett syndrome	<i>MECP2</i>	Methyl CpG binding protein 2	XL	88% females with Rett syndrome <sup>13</sup>
	<i>CDKL5</i>	Cyclin-dependent kinase-like 5	XL	2-8% females with atypical Rett syndrome <sup>14,15</sup>
	<i>CTNNB1</i>	Catenin beta 1	AD	Rare in Rett-like syndromes <sup>29</sup>

	<i>DDX3X</i>	DEAD-box helicase 3, X-linked	XL	Rare in Rett-like syndromes <sup>30</sup>
	<i>FOXG1*</i>	Forkhead box protein G1	AD	~1% Rett syndrome overall <sup>16</sup> ; 25% congenital variant of Rett <sup>17</sup>
	<i>GABBR2</i>	Gamma-aminobutyric acid type B receptor subunit 2	AD	Rare in Rett-like syndromes <sup>31</sup>
	<i>IQSEC2</i>	IQ motif and sec7 domain 7	AD	Rare in Rett-like syndromes <sup>39</sup>
	<i>KCNA2</i>	Potassium voltage-gated channel subfamily A member 2	AD	Rare in Rett-like syndromes <sup>32</sup>
	<i>MBD5</i>	Methyl-CpG-binding domain protein 5	AD	Rare in Rett-like syndromes <sup>26</sup>
	<i>MEF2C</i>	Myocyte-specific enhancer factor 2C	AD	Rare in Rett-like syndromes <sup>23</sup>
	<i>PCDH19**</i>	Protocadherin-19	XL	Rare in Rett-like syndromes <sup>40</sup>
	<i>SATB2</i>	SATB homeobox 2	AD	Rare in Rett-like syndromes <sup>33</sup>
	<i>SHANK3**</i>	SH3 and multiple ankyrin repeat domains protein 3	AD	Rare in Rett-like syndromes <sup>41</sup>
	<i>STXBP1</i>	Syntaxin binding protein 1	AD	Rare in Rett-like syndromes <sup>34</sup>
	<i>TBL1XR1</i>	Transducin (beta)-like 1 X-linked receptor 1	AD	Rare in Rett-like syndromes <sup>35</sup>
Angelman/ Angelman-like syndromes	<i>WDR45</i>	WD repeat domain 45	XL	Rare in Rett-like syndromes <sup>27,28</sup>
	<i>UBE3A***</i>	Ubiquitin protein ligase E3A	AD	78% Angelman syndrome by methylation <sup>18</sup> ; 11% Angelman syndrome by sequencing of UBE3A <sup>18</sup>
	<i>SLC9A6</i>	Sodium/hydrogen exchanger 6	XL	6% Angelman-like syndrome <sup>7</sup>
	<i>TCF4</i>	Transcription factor 4	AD	2% Angelman syndrome <sup>19</sup>
	<i>DYRK1A</i>	Dual-specificity tyrosine-(Y)-phosphorylation regulated kinase 1A	AD	Rare <sup>37</sup>
	<i>EHMT1</i>	Histone-lysine N-methyltransferase EHMT1	AD	~75% of patients with Kleefstra syndrome, including microdeletion of 9q34.3 <sup>38</sup>
Pitt-Hopkins syndrome (PHS)	<i>MBD5</i>	Methyl-CpG-binding domain protein 5	AD	Rare in Angelman-like syndrome <sup>25</sup>
	<i>TCF4</i>	Transcription factor 4	AD	36% PHS <sup>19</sup>
	<i>NRXN1</i>	Neurexin-1	AR	Rare in PHS <sup>20</sup>
Mowat-Wilson syndrome	<i>CNTNAP2</i>	Contactin-associated protein-like 2	AR	Rare in PHS <sup>20</sup>
	<i>ZEB2</i>	Zinc finger E-box-binding homeobox 2	AD	95% Mowat Wilson syndrome <sup>21,22</sup>

Alpha-thalassemia X-linked intellectual disability	ATRX	Transcriptional regulator ATRX	XL	95% males with alpha-thalassemia X-linked intellectual disability <sup>36</sup>
--	------	--------------------------------	----	---

\* This panel does not include deletion/duplication testing of FOXP1

\*\* Only whole gene deletion/duplications may be detected for the PCDH19 and SHANK3 genes

\*\*\* This panel includes MS-MLPA to evaluate for abnormal methylation of the UBE3A gene

## REFERENCES:

1. Christodoulou J and Ho G. MECP2- Related Disorders. In: GeneReviews at GeneTests: Medical Genetics Information Resource (database online). Copyright, University of Washington, Seattle. 1997-2011. Available at <http://www.genetests.org>. Accessed February 2012.
2. Mencarelli et al., (2010) J Med Genet 47:49-53.
3. Le Guen et al., (2011) Neurogenet 12:1-8.
4. Williams, C.A. and Dagli A.I. (Updated June 2011). Angelman Syndrome. In: GeneReviews at GeneTests: Medical Genetics Information Resource (database online). Copyright, University of Washington, Seattle. 1997-2009. Available at <http://www.genetests.org>. Accessed July 2012.
5. Van Buggenhout and Frynes (2009) Eur J Hum Genet 17:1367-1373.
6. Schroer et al., (2009) Am J Med Genet 152A:2775-2783.
7. Gillfillan et al., (2008) Am J Hum Genet 82(4):1003-1010.
8. Christianson et al., (1999) Am J. Hum Genet. 36:759-766.
9. Zweier et al., (2008) J Med Genet 45:738-744.
10. Marangi et al., (2011) Am J Med Genet A 155A:1536-1545.
11. Whalen et al., (2012) Hum Mutat 33:64-72.
12. Garavelli et al., (2009) Am J Med Genet A 149A:417-426.
13. Li et al., (2006) J Hum Genet 52:38-47.
14. Tao et al., (2004) Am J Hum Genet 75:1149-1154.
15. Rosas-Vargas et al., (2008) J Med Genet 45:172-178.
16. Bahi-Buisson et al., (2010) 11:241-249.
17. Mencarelli et al., (2010) J Med Genet 47:49-53.
18. Lossie, et al, (2001) J Med Genet 38(12):834-845.
19. de Pontual et al., (2009) Hum Mutat 30:669-676.
20. Zweier et al., (2009) Am J Hum Genet 85:655-666.
21. Wilson et al., (2003) Am J Med Genet A 119A:257-265.
22. Zweier et al., (2005) Eur J Med Genet 48:97-111.
23. Zweier et al., (2010) Hum Mutat 31:722-733.
24. Bienvenu et al., (2013) Neurogenet 14:71-75.
25. Hodge et al., (2013) Mol Psych Apr 16. doi: 10.1038/mp.2013.42.
26. Jaillard et al., (2009) J Med Genet 46:847-855.
27. Depienne et al., (2009) Plos Genet 5(2):e1000381.
28. Marini et al., (2010) Neurology 75:646-653.
29. Percy et al. (2018) Transl Sci Rare Dis 3 (1):49-53 (PMID: 29682453).
30. Snijders et al. (2015) Am. J. Hum. Genet. 97 (2):343-52 (PMID: 26235985).
31. Yoo et al. (2017) Ann. Neurol. 82 (3):466-478 (PMID: 28856709).
32. Allou et al. (2017) Clin. Genet. 91 (3):431-440 (PMID: 27062609).
33. Zarate et al. (2018) Am. J. Med. Genet. A 176 (4):925-935 (PMID: 29436146)
34. Romaniello et al. (2015) Neuroreport 26 (5):254-7 (PMID: 25714420).
35. Saitsu et al. (2014) Journal Of Human Genetics 59 (10):581-3 (PMID: 25102098)
36. Stevenson RE. Alpha-Thalassemia X-Linked Intellectual Disability Syndrome. 2000 Jun 19 [Updated 2014 Nov 6]. In: Adam MP, Ardinger HH, Pagon RA, et al., editors. GeneReviews® [Internet]. Seattle (WA): University of Washington, Seattle; 1993-2018. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK1449/>.
37. Courcet et al. (2012) Journal Of Medical Genetics 49 (12):731-6 (PMID: 23099646)
38. Willemsen et al. (2012) Mol Syndromol 2 (3-5):202-212 (PMID: 22670141)
39. Tran et al. (2014) Eur. J. Hum. Genet. 22 (2):289-92 (PMID: 23674175)
40. Depienne et al., (2010) Hum Mutat 32:E1959-E1975
41. Hara et al. (2015) Am. J. Med. Genet. A 167 (7):1593-6 (PMID: 25931020)