

## OncoGeneDx: Comprehensive Common Cancer Panel

### PANEL GENE LIST

*APC, ATM, AXIN2, BAP1, BARD1, BMPR1A, BRCA1, BRCA2, BRIP1, CDH1, CDK4, CDKN2A, CHEK2, CTNNA1, EPCAM\*, FANCC, FANCM, FH, FLCN, HOXB13, MET, MITF, MLH1, MSH2, MSH3, MSH6, MUTYH, NBN, NF1, NTHL1, PALB2, PMS2, POLD1, POLE, POT1, PTEN, RAD51C, RAD51D, RECQL, SCG5/GREM1\*, SDHB, SDHC, SDHD, SMAD4, STK11, TP53, TSC1, TSC2, VHL*

\*Testing includes sequencing and deletion/duplication analysis for all genes except *EPCAM* and *SCG5/GREM1* (del/dup only).

### CLINICAL FEATURES

Cancer is a common disease affecting approximately 1 in 3 individuals in the U.S.<sup>1</sup> While the majority of cancers are sporadic in nature, some families have hereditary forms of cancer that are associated with increased cancer risks compared with the general population.

Approximately 5-10% of cancer cases are thought to be due to a hereditary predisposition. The features of a personal and/or family history of cancer that are suggestive of a hereditary cancer predisposition include: young ages at diagnosis, multiple primary cancers in a single individual, and several relatives affected with the same type of cancer or related cancers spanning multiple generations.

For some of the well-described hereditary conditions discussed below, clinical diagnostic criteria based on personal medical history and family history are available to help identify patients most likely to have a hereditary cancer syndrome. In many cases, however, patients do not meet the clinical diagnostic criteria or the criteria may overlap for multiple conditions, making it difficult to decide which genes should be tested and in what order. The OncoGeneDx Comprehensive Common Cancer Panel offered at GeneDx includes analysis of 49 genes associated with hereditary predisposition to various cancers including the most common hereditary cancer syndromes such as Hereditary Breast and Ovarian Cancer syndrome (*BRCA1, BRCA2*) and Lynch syndrome (*MLH1, MSH2, MSH6, PMS2, EPCAM*), as well as newly described genes.

Many of the genes on this panel are involved in the mismatch repair pathway, the Fanconi anemia pathway and/or play a role in DNA damage repair. Newer genes that have been identified in families with cancer have been included in the panel to make it as comprehensive as possible. These genes include, but are not limited to, *AXIN2, CTNNA1, MSH3, NTHL1, and RECQL*. The evidence available to date may be derived from a small number of patients with wide confidence intervals or is based upon an ethnic cohort with one specific variant. Accurate risk assessment may be complicated by the low penetrance of pathogenic variants in these genes and/or ascertainment bias.

### INHERITANCE PATTERN

Most genes on this panel are associated with an autosomal dominant cancer risk with the exception of *MSH3, MUTYH* and *NTHL1*, which are associated with an autosomal recessive cancer risk. Some of the genes on this panel are also associated with extremely rare conditions when inherited in an autosomal recessive fashion. The specifics of this inheritance are outlined in the table below.

### TEST METHODS

Genomic DNA is extracted from the submitted specimen. For skin punch biopsies, fibroblasts are cultured and used for DNA extraction. The DNA is enriched for the complete coding regions and splice junctions of the genes on this panel using a proprietary targeted capture system developed by GeneDx for next-generation sequencing with CNV calling (NGS-CNV). For *PTEN* nucleotides c.-700 through c.-1300 in the promoter region, and for *APC*, promoters 1A and 1B are also captured. The enriched targets are simultaneously sequenced with paired end reads on an Illumina platform. Bi-directional sequence reads are assembled and aligned to reference sequences

based on NCBI RefSeq transcripts and human genome build GRCh37/UCSC hg19. After gene specific filtering, data are analyzed to identify sequence variants and most deletions and duplications involving coding exons; however, technical limitations and inherent sequence properties effectively reduce this resolution for some genes. Concurrent *MSH2* Exons 1-7 Inversion analysis from NGS data is also performed. For *EPCAM* and *SCG5*, deletion/duplication analysis, but not sequencing, is performed. Alternative sequencing or copy number detection methods are used to analyze or confirm regions with inadequate sequence or copy number data by NGS. Reportable variants include pathogenic variants, likely pathogenic variants and variants of uncertain significance. Likely benign and benign variants, if present, are not routinely reported but are available upon request.

## TEST SENSITIVITY

The clinical sensitivity of sequencing and deletion/duplication analysis of the 49 genes included in the OncoGeneDx Comprehensive Common Cancer Panel depends in part on the patient's clinical phenotype and family history. In general, the sensitivity is highest for individuals with features suggestive of a hereditary predisposition to cancer as outlined above. DNA sequencing will detect nucleotide substitutions and small insertions and deletions, while NGS CNV analysis, array CGH, or MLPA will detect exon-level deletions and duplications. These methods are expected to be greater than 99% sensitive in detecting pathogenic variants identifiable by sequencing or CNV technology.

Genetic testing using the methods applied at GeneDx is expected to be highly accurate. Normal findings do not rule out the diagnosis of a genetic disorder since some genetic abnormalities may be undetectable by this test. The methods used cannot reliably detect deletions of 20bp to 250bp in size, or insertions of 10bp to 250 bp in size. Sequencing cannot detect low-level mosaicism. The copy number assessment methods used with this test cannot reliably detect mosaicism and cannot identify balanced chromosome aberrations. Rarely, incidental findings of large chromosomal rearrangements outside the gene of interest may be identified. Regions of certain genes have inherent sequence properties (for example: repeat, homology, or pseudogene regions, high GC content, rare polymorphisms) that yield suboptimal data, potentially impairing accuracy of the results. False negatives may also occur in the setting of bone marrow transplantation, recent blood transfusion, or suboptimal DNA quality. In individuals with active or chronic hematologic neoplasms or conditions, there is a possibility that testing may detect an acquired somatic variant, resulting in a false positive result. As the ability to detect genetic variants and naming conventions can differ among laboratories, rare false negative results may occur when no positive control is provided for testing of a specific variant identified at another laboratory. The chance of a false positive or false negative result due to laboratory errors incurred during any phase of testing cannot be completely excluded. Interpretations are made with the assumption that any clinical information provided, including family relationships, are accurate. Consultation with a genetics professional is recommended for interpretation of results.

Gene	Protein	Inheritance	Disease Associations
<i>APC</i> <sup>2,3</sup>	ADENOMATOUS POLYPOSIS COLI PROTEIN	AD	Familial adenomatous polyposis (FAP)-associated condition: colorectal, duodenal or periampullary, gastric, thyroid, pancreatic, brain (medulloblastoma) & liver (hepatoblastoma) cancers, desmoid tumors, gastrointestinal polyps
<i>ATM</i> <sup>4-8</sup>	SERINE-PROTEIN KINASE ATM	AD	Breast, pancreatic, prostate & colon cancer
		AR	Ataxia telangiectasia

Gene	Protein	Inheritance	Disease Associations
<i>AXIN2</i> <sup>9,10</sup>	AXIN-2	AD	Colon cancer, colon polyps
<i>BAP1</i> <sup>11,12</sup>	UBIQUITIN CARBOXYL-TERMINAL HYDROLASE BAP1	AD	Uveal/cutaneous melanoma, basal cell carcinoma, mesothelioma, renal cancer
<i>BARD1</i> <sup>13-15</sup>	BRCA1-ASSOCIATED RING DOMAIN PROTEIN 1	AD	Breast cancer
<i>BMPR1A</i> <sup>16,17</sup>	BONE MORPHOGENETIC PROTEIN RECEPTOR TYPE-1A	AD	Juvenile Polyposis syndrome (JPS): colorectal, gastric (if gastric polyps), small bowel & pancreatic cancer, gastrointestinal polyps
<i>BRCA1</i> <sup>18,19</sup>	BREAST CANCER TYPE 1 SUSCEPTIBILITY PROTEIN	AD	Hereditary breast and ovarian cancer (HBOC) syndrome: breast, ovarian, pancreatic & prostate cancer
<i>BRCA2</i> <sup>18,20</sup>	BREAST CANCER TYPE 2 SUSCEPTIBILITY PROTEIN	AD	Hereditary breast and ovarian cancer (HBOC) syndrome: breast, ovarian, pancreatic, prostate cancer & melanoma
		AR	Fanconi anemia
<i>BRIP1</i> <sup>13,20-22</sup>	FANCONI ANEMIA GROUP J PROTEIN	AD	Breast, ovarian & prostate cancer
		AR	Fanconi anemia
<i>CDH1</i> <sup>23-25</sup>	CADHERIN 1	AD	Hereditary diffuse gastric cancer (HDGC) syndrome: gastric (diffuse), breast & colon (signet ring) cancer
<i>CDK4</i> <sup>26-28</sup>	CYCLIN-DEPENDENT KINASE 4	AD	Melanoma, non-melanoma skin & pancreatic cancer
<i>CDKN2A</i> <sup>28-33</sup>	CYCLIN-DEPENDENT KINASE INHIBITOR 2A, TUMOR SUPPRESSOR ARF	AD	Familial atypical multiple mole melanoma (FAMMM) syndrome: melanoma, pancreatic cancer & astrocytoma
<i>CHEK2</i> <sup>15,34-37</sup>	SERINE/THREONINE-PROTEIN KINASE CHK2	AD	Breast, colon, prostate, gastric & thyroid cancer
<i>CTNNA1</i> <sup>25,38</sup>	CATENIN ALPHA-1	AD	Diffuse gastric cancer
<i>EPCAM</i> <sup>39-41</sup>	EPITHELIAL CELL ADHESION MOLECULE	AD	Lynch syndrome (LS): colorectal, endometrial, ovarian, gastric, pancreatic, biliary tract, urinary tract, small bowel, prostate & brain cancer, sebaceous neoplasms
		AR	Constitutional mismatch repair deficiency syndrome
<i>FANCC</i> <sup>20,42,43</sup>	FANCONI ANEMIA GROUP C PROTEIN	AD	Breast cancer
		AR	Fanconi anemia
<i>FANCM</i> <sup>15,44</sup>	FANCONI ANEMIA GROUP M PROTEIN	AD	Breast cancer
		AR	Fanconi anemia-like cancer susceptibility
<i>FH</i> <sup>45,46</sup>	FUMARATE HYDRATASE, MITOCHONDRIAL	AD	Hereditary leiomyomatosis and renal cell cancer (HLRCC): renal cancer, leiomyomas, pheochromocytoma, paraganglioma
		AR	Fumarate hydratase deficiency

Gene	Protein	Inheritance	Disease Associations
<i>FLCN</i> <sup>47</sup>	FOLLICULIN	AD	Birt-Hogg-Dubé syndrome (BHD): renal cancer
<i>HOXB13</i> <sup>48,49</sup>	HOMEBOX PROTEIN HOX-B13	AD	Prostate cancer
<i>MET</i> <sup>50,51</sup>	HEPATOCTE GROWTH FACTOR RECEPTOR	AD	Hereditary papillary renal carcinoma (HPRC): renal cancer (type I papillary)
<i>MITF</i> <sup>52,53</sup>	MICROPHthalmia-ASSOCIATED TRANSCRIPTION FACTOR	AD	Renal cancer, melanoma
<i>MLH1</i> <sup>54,55</sup>	DNA MISMATCH REPAIR PROTEIN MLH1	AD	Lynch syndrome (LS): colorectal, endometrial, ovarian, gastric, pancreatic, biliary tract, urinary tract, small bowel, prostate & brain cancer, sebaceous neoplasms
		AR	Constitutional mismatch repair deficiency syndrome
<i>MSH2</i> <sup>54,55</sup>	DNA MISMATCH REPAIR PROTEIN MSH2	AD	Lynch syndrome (LS): colorectal, endometrial, ovarian, gastric, pancreatic, biliary tract, urinary tract, small bowel, prostate & brain cancer, sebaceous neoplasms
		AR	Constitutional mismatch repair deficiency syndrome
<i>MSH3</i> <sup>56,57</sup>	DNA MISMATCH REPAIR PROTEIN MSH3	AR	Colorectal cancer, colonic polyposis
<i>MSH6</i> <sup>54,55</sup>	DNA MISMATCH REPAIR PROTEIN MSH6	AD	Lynch syndrome (LS): colorectal, endometrial, ovarian, gastric, pancreatic, biliary tract, urinary tract, small bowel, prostate & brain cancer, sebaceous neoplasms
		AR	Constitutional mismatch repair deficiency syndrome
<i>MUTYH</i> <sup>58,59</sup>	ADENINE DNA GLYCOSYLASE	AR	<i>MUTYH</i> -associated polyposis (MAP): colorectal, small bowel & endometrial serous cancer, gastrointestinal polyps
<i>NBN</i> <sup>60-62</sup>	NIBRIN	AD	Prostate cancer
		AR	Nijmegen breakage syndrome
<i>NF1</i> <sup>63,64</sup>	NEUROFIBROMIN	AD	Neurofibromatosis type 1 (NF1) syndrome: breast cancer, GIST, optic nerve gliomas, pheochromocytoma, MPNST, neurofibromas, brain tumors
<i>NTHL1</i> <sup>57,65</sup>	ENDONUCLEASE III-LIKE 1	AR	Colon cancer, colon polyps
<i>PALB2</i> <sup>13,20,66,67</sup>	PARTNER AND LOCALIZER OF BRCA2	AD	Breast, pancreatic, ovarian & prostate cancer
		AR	Fanconi anemia
<i>PMS2</i> <sup>54,55</sup>	MISMATCH REPAIR ENDONUCLEASE PMS2	AD	Lynch syndrome (LS): colorectal, endometrial, ovarian, gastric, pancreatic,

Gene	Protein	Inheritance	Disease Associations
			biliary tract, urinary tract, small bowel, prostate & brain cancer, sebaceous neoplasms
		AR	Constitutional mismatch repair deficiency syndrome
<i>POLD1</i> <sup>68,69</sup>	DNA POLYMERASE DELTA CATALYTIC SUBUNIT	AD	Colon & endometrial cancer, colon polyps
		AD	Colon cancer, gastrointestinal polyps
<i>POLE</i> <sup>70,71</sup>	DNA POLYMERASE EPSILON CATALYTIC SUBUNIT A	AR	Intrauterine growth restriction, metaphyseal dysplasia, congenital adrenal hypoplasia, and genitourinary anomalies in males (IMAGe) with variable immunodeficiency
<i>POT1</i> <sup>72,73</sup>	PROTECTION OF TELOMERES 1	AD	Melanoma, hematologic malignancy, brain glial tumors
<i>PTEN</i> <sup>74</sup>	PHOSPHATIDYLINOSITOL 3,4,5-TRISPHOSPHATE 3-PHOSPHATASE AND DUAL-SPECIFICITY PROTEIN PHOSPHATASE PTEN	AD	PTEN hamartoma tumor syndrome (PHTS): breast, thyroid, endometrial, colon, melanoma, renal cancer, gastrointestinal polyps, Lhermitte-Duclos Disease
<i>RAD51C</i> <sup>15,21,75</sup>	DNA REPAIR PROTEIN RAD51 HOMOLOG 3	AD	Breast, ovarian & prostate cancer
		AR	Fanconi anemia
<i>RAD51D</i> <sup>15,21,75</sup>	DNA REPAIR PROTEIN RAD51 HOMOLOG 4	AD	Breast, ovarian & prostate cancer
<i>RECQL</i> <sup>15,76</sup>	RECQ PROTEIN-LIKE	AD	Breast cancer
<i>SCG5/ GREM1</i> <sup>77-79</sup>	NEUROENDOCRINE PROTEIN 7B2/GREMLIN-1	AD	Hereditary mixed polyposis syndrome (HMPS): colon cancer, colon polyps
<i>SDHB</i> <sup>80,81</sup>	SUCCINATE DEHYDROGENASE [UBIQUINONE] IRON-SULFUR SUBUNIT, MITOCHONDRIAL	AD	Hereditary paraganglioma/ pheochromocytoma (PGL/PCC) syndrome: paraganglioma, pheochromocytoma, renal cancer, GIST
		AR	Isolated complex II deficiency
<i>SDHC</i> <sup>80</sup>	SUCCINATE DEHYDROGENASE CYTOCHROME B560 SUBUNIT, MITOCHONDRIAL	AD	Hereditary paraganglioma/ pheochromocytoma (PGL/PCC) syndrome: paraganglioma, pheochromocytoma, renal cancer, GIST
<i>SDHD</i> <sup>80,82</sup>	SUCCINATE DEHYDROGENASE [UBIQUINONE] CYTOCHROME B SMALL SUBUNIT, MITOCHONDRIAL	AD	Hereditary paraganglioma/ pheochromocytoma (PGL/PCC) syndrome: paraganglioma, pheochromocytoma, renal cancer, GIST
		AR	Isolated complex II deficiency
<i>SMAD4</i> <sup>83,84</sup>	MOTHERS AGAINST DECAPENTAPLEGIC HOMOLOG 4	AD	Juvenile polyposis syndrome (JPS): colorectal, gastric (if gastric polyps), small bowel & pancreatic cancer, gastrointestinal polyps

Gene	Protein	Inheritance	Disease Associations
<i>STK11</i> <sup>85-87</sup>	SERINE/THREONINE-PROTEIN KINASE STK11	AD	Peutz-Jeghers syndrome (PJS): breast, colorectal, pancreatic, gastric, small bowel, ovarian, lung, cervical & endometrial cancer, testicular tumors (LCCSCT), gastrointestinal polyps
<i>TP53</i> <sup>88</sup>	CELLULAR TUMOR ANTIGEN P53	AD	Li-Fraumeni syndrome (LFS): breast cancer, sarcoma, brain cancer, hematologic malignancies, adrenocortical carcinoma, among others**
<i>TSC1</i> <sup>89-91</sup>	HAMARTIN	AD	Tuberous sclerosis complex (TSC): renal cancer/tumors, central nervous system tumors, cardiac rhabdomyomas, angiomyolipomas
<i>TSC2</i> <sup>89-91</sup>	TUBERIN	AD	Tuberous sclerosis complex (TSC): renal cancer/tumors, central nervous system tumors, cardiac rhabdomyomas, angiomyolipomas
<i>VHL</i> <sup>92,93</sup>	VON HIPPEL-LINDAU DISEASE TUMOR SUPPRESSOR	AD	von Hippel-Lindau (VHL) disease: renal cancer (clear cell), pancreatic neuroendocrine tumors, hemangioblastoma, pheochromocytoma, endolymphatic sac tumors

Because of evolving and expanding phenotypes, this list of cancer/tumor types is not exhaustive. Gene-specific risk for some of the cancers and other features listed are not well-defined.

\*\* High overall risk of cancer: 75% lifetime risk for males to develop cancer, nearly 100% risk for females.

#### Abbreviations:

AD – Autosomal Dominant

AR – Autosomal Recessive

CGH – Comparative genomic hybridization

GIST – Gastrointestinal stromal tumor

LCCSCT - Large cell-calcifying Sertoli cell tumors

MLPA – Multiplex ligation-dependent probe amplification

MPNST - Malignant peripheral nerve sheath tumors

#### REFERENCES:

1. Probability of Developing or Dying of Cancer - Surveillance Research Program. <https://surveillance.cancer.gov/devcan/>.
2. Jasperson, K. W., Patel, S. G. & Ahnen, D. J. APC-Associated Polyposis Conditions. in *GeneReviews*® (eds. Adam, M. P. et al.) (University of Washington, Seattle, 1993).
3. Li, J. *et al.* Point Mutations in Exon 1B of APC Reveal Gastric Adenocarcinoma and Proximal Polyposis of the Stomach as a Familial Adenomatous Polyposis Variant. *Am J Hum Genet* **98**, 830–842 (2016).
4. Gatti, R. & Perlman, S. Ataxia-Telangiectasia. in *GeneReviews*® (eds. Adam, M. P. et al.) (University of Washington, Seattle, 1993).
5. Hu, C. *et al.* Association Between Inherited Germline Mutations in Cancer Predisposition Genes and Risk of Pancreatic Cancer. *JAMA* **319**, 2401–2409 (2018).
6. Lu, H.-M. *et al.* Association of Breast and Ovarian Cancers With Predisposition Genes Identified by Large-Scale Sequencing. *JAMA Oncol* **5**, 51–57 (2019).
7. LaDuca, H. *et al.* A clinical guide to hereditary cancer panel testing: evaluation of gene-specific cancer associations and sensitivity of genetic testing criteria in a cohort of 165,000 high-risk patients. *Genet Med* **22**, 407–415 (2020).
8. Hall, M. J. *et al.* Germline Pathogenic Variants in the Ataxia Telangiectasia Mutated (ATM) Gene are Associated with High and Moderate Risks for Multiple Cancers. *Cancer Prev Res (Phila)* **14**, 433–440 (2021).
9. Marvin, M. L. *et al.* AXIN2-associated autosomal dominant ectodermal dysplasia and neoplastic syndrome. *Am. J. Med. Genet. A* **155A**, 898–902 (2011).
10. Lammi, L. *et al.* Mutations in AXIN2 cause familial tooth agenesis and predispose to colorectal cancer. *Am. J. Hum. Genet.* **74**, 1043–1050 (2004).
11. Pilarski, R., Rai, K., Cebulla, C. & Abdel-Rahman, M. BAP1 Tumor Predisposition Syndrome. in *GeneReviews*® (eds. Adam, M. P. et al.) (University of Washington, Seattle, 1993).
12. Rai, K., Pilarski, R., Cebulla, C. M. & Abdel-Rahman, M. H. Comprehensive review of BAP1 tumor predisposition syndrome with report of two new cases. *Clin. Genet.* **89**, 285–294 (2016).
13. Shimelis, H. *et al.* Triple-Negative Breast Cancer Risk Genes Identified by Multigene Hereditary Cancer Panel Testing. *J Natl Cancer Inst* **110**, 855–862 (2018).

14. Suszynska, M. & Kozlowski, P. Summary of BARD1 Mutations and Precise Estimation of Breast and Ovarian Cancer Risks Associated with the Mutations. *Genes (Base)* **11**, E798 (2020).
15. Breast Cancer Association Consortium *et al.* Breast Cancer Risk Genes - Association Analysis in More than 113,000 Women. *N Engl J Med* **384**, 428–439 (2021).
16. Syngal, S. *et al.* ACG clinical guideline: Genetic testing and management of hereditary gastrointestinal cancer syndromes. *Am J Gastroenterol* **110**, 223–262; quiz 263 (2015).
17. Brosens, L. A. A. *et al.* Risk of colorectal cancer in juvenile polyposis. *Gut* **56**, 965–967 (2007).
18. Petrucelli, N., Daly, M. B. & Pal, T. BRCA1- and BRCA2-Associated Hereditary Breast and Ovarian Cancer. in *GeneReviews®* (eds. Adam, M. P. *et al.*) (University of Washington, Seattle, 1993).
19. Freire, B. L. *et al.* Homozygous loss of function BRCA1 variant causing a Fanconi-anemia-like phenotype, a clinical report and review of previous patients. *Eur J Med Genet* **61**, 130–133 (2018).
20. Mehta, P. A. & Ebens, C. Fanconi Anemia. in *GeneReviews®* (eds. Adam, M. P. *et al.*) (University of Washington, Seattle, 1993).
21. Pritchard, C. C. *et al.* Inherited DNA-Repair Gene Mutations in Men with Metastatic Prostate Cancer. *N. Engl. J. Med.* **375**, 443–453 (2016).
22. Rafnar, T. *et al.* Mutations in BRIP1 confer high risk of ovarian cancer. *Nat Genet* **43**, 1104–1107 (2011).
23. Xicola, R. M. *et al.* Clinical features and cancer risk in families with pathogenic CDH1 variants irrespective of clinical criteria. *J Med Genet* **56**, 838–843 (2019).
24. Roberts, M. E. *et al.* Comparison of CDH1 Penetrance Estimates in Clinically Ascertained Families vs Families Ascertained for Multiple Gastric Cancers. *JAMA Oncol* **5**, 1325–1331 (2019).
25. Blair, V. R. *et al.* Hereditary diffuse gastric cancer: updated clinical practice guidelines. *Lancet Oncol* **21**, e386–e397 (2020).
26. Goldstein, A. M. *et al.* High-risk melanoma susceptibility genes and pancreatic cancer, neural system tumors, and uveal melanoma across GenoMEL. *Cancer Res.* **66**, 9818–9828 (2006).
27. Puntervoll, H. E. *et al.* Melanoma prone families with CDK4 germline mutation: phenotypic profile and associations with MC1R variants. *J. Med. Genet.* **50**, 264–270 (2013).
28. Leachman, S. A. *et al.* Identification, genetic testing, and management of hereditary melanoma. *Cancer Metastasis Rev.* **36**, 77–90 (2017).
29. Goldstein, A. M. *et al.* Features associated with germline CDKN2A mutations: a GenoMEL study of melanoma-prone families from three continents. *J. Med. Genet.* **44**, 99–106 (2007).
30. Kefford, R. F., Newton Bishop, J. A., Bergman, W. & Tucker, M. A. Counseling and DNA testing for individuals perceived to be genetically predisposed to melanoma: A consensus statement of the Melanoma Genetics Consortium. *J. Clin. Oncol.* **17**, 3245–3251 (1999).
31. Czajkowski, R., Placek, W., Drewa, G., Czajkowska, A. & Uchańska, G. FAMMM syndrome: pathogenesis and management. *Dermatol Surg* **30**, 291–296 (2004).
32. Canto, M. I. *et al.* International Cancer of the Pancreas Screening (CAPS) Consortium summit on the management of patients with increased risk for familial pancreatic cancer. *Gut* **62**, 339–347 (2013).
33. Soua, E., Eliades, P. J., Shannon, K., Stratigos, A. J. & Tsao, H. Hereditary melanoma: Update on syndromes and management: Emerging melanoma cancer complexes and genetic counseling. *J. Am. Acad. Dermatol.* **74**, 411–420; quiz 421–422 (2016).
34. Hu, C. *et al.* A Population-Based Study of Genes Previously Implicated in Breast Cancer. *N Engl J Med* **384**, 440–451 (2021).
35. Couch, F. J. *et al.* Associations Between Cancer Predisposition Testing Panel Genes and Breast Cancer. *JAMA Oncol* **3**, 1190–1196 (2017).
36. Siolek, M. *et al.* CHEK2 mutations and the risk of papillary thyroid cancer. *Int J Cancer* **137**, 548–552 (2015).
37. Han, F., Guo, C. & Liu, L. The effect of CHEK2 variant I157T on cancer susceptibility: evidence from a meta-analysis. *DNA Cell Biol.* **32**, 329–335 (2013).
38. Clark, D. F. *et al.* Loss-of-function variants in CTNNA1 detected on multigene panel testing in individuals with gastric or breast cancer. *Genet Med* **22**, 840–846 (2020).
39. Pathak, S. J. *et al.* EPCAM mutation update: Variants associated with congenital tufting enteropathy and Lynch syndrome. *Hum Mutat* **40**, 142–161 (2019).
40. Cini, G. *et al.* Toward a better definition of EPCAM deletions in Lynch Syndrome: Report of new variants in Italy and the associated molecular phenotype. *Mol Genet Genomic Med* **7**, e587 (2019).
41. Kempers, M. J. E. *et al.* Risk of colorectal and endometrial cancers in EPCAM deletion-positive Lynch syndrome: a cohort study. *Lancet Oncol* **12**, 49–55 (2011).
42. Thompson, E. R. *et al.* Exome sequencing identifies rare deleterious mutations in DNA repair genes FANCC and BLM as potential breast cancer susceptibility alleles. *PLoS Genet.* **8**, e1002894 (2012).
43. Berwick, M. *et al.* Genetic heterogeneity among Fanconi anemia heterozygotes and risk of cancer. *Cancer Res.* **67**, 9591–9596 (2007).
44. Catucci, I. *et al.* Individuals with FANCM biallelic mutations do not develop Fanconi anemia, but show risk for breast cancer, chemotherapy toxicity and may display chromosome fragility. *Genet Med* **20**, 452–457 (2018).
45. Kamihara, J., Schultz, K. A. & Rana, H. Q. FH Tumor Predisposition Syndrome. in *GeneReviews®* (eds. Adam, M. P. *et al.*) (University of Washington, Seattle, 1993).
46. Coman, D., Kranc, K. R. & Christodoulou, J. Fumarate Hydratase Deficiency. in *GeneReviews®* (eds. Adam, M. P. *et al.*) (University of Washington, Seattle, 1993).
47. Menko, F. H. *et al.* Birt-Hogg-Dubé syndrome: diagnosis and management. *Lancet Oncol.* **10**, 1199–1206 (2009).
48. Karlsson, R. *et al.* A population-based assessment of germline HOXB13 G84E mutation and prostate cancer risk. *Eur Urol* **65**, 169–176 (2014).
49. Nyberg, T. *et al.* Homeobox B13 G84E Mutation and Prostate Cancer Risk. *Eur Urol* **75**, 834–845 (2019).
50. Lubensky, I. A. *et al.* Hereditary and sporadic papillary renal carcinomas with c-met mutations share a distinct morphological phenotype. *Am. J. Pathol.* **155**, 517–526 (1999).
51. Schmidt, L. S. *et al.* Early onset hereditary papillary renal carcinoma: germline missense mutations in the tyrosine kinase domain of the met proto-oncogene. *J. Urol.* **172**, 1256–1261 (2004).
52. Bertolotto, C. *et al.* A SUMOylation-defective MITF germline mutation predisposes to melanoma and renal carcinoma. *Nature* **480**, 94–98 (2011).
53. Ghiorzo, P. *et al.* Prevalence of the E318K MITF germline mutation in Italian melanoma patients: associations with histological subtypes and family cancer history. *Pigment Cell Melanoma Res* **26**, 259–262 (2013).
54. Idos, G. & Valle, L. Lynch Syndrome. in *GeneReviews®* (eds. Adam, M. P. *et al.*) (University of Washington, Seattle, 1993).
55. Wimmer, K. *et al.* Diagnostic criteria for constitutional mismatch repair deficiency syndrome: suggestions of the European consortium ‘care for CMMRD’ (C4CMMRD). *J. Med. Genet.* **51**, 355–365 (2014).

56. Adam, R. *et al.* Exome Sequencing Identifies Biallelic MSH3 Germline Mutations as a Recessive Subtype of Colorectal Adenomatous Polyposis. *Am J Hum Genet* **99**, 337–351 (2016).
57. Terradas, M. *et al.* Contribution to colonic polyposis of recently proposed predisposing genes and assessment of the prevalence of NTHL1- and MSH3-associated polyposes. *Hum Mutat* **40**, 1910–1923 (2019).
58. Win, A. K. *et al.* Risk of extracolonic cancers for people with biallelic and monoallelic mutations in MUTYH. *Int J Cancer* **139**, 1557–1563 (2016).
59. Win, A. K. *et al.* Risk of colorectal cancer for carriers of mutations in MUTYH, with and without a family history of cancer. *Gastroenterology* **146**, 1208–1211.e1–5 (2014).
60. Varon, R., Demuth, I. & Chrzanowska, K. H. Nijmegen Breakage Syndrome. in *GeneReviews®* (eds. Adam, M. P. et al.) (University of Washington, Seattle, 1993).
61. Wokolorczyk, D. *et al.* Mutations in ATM, NBN and BRCA2 predispose to aggressive prostate cancer in Poland. *Int J Cancer* **147**, 2793–2800 (2020).
62. Cybulski, C. *et al.* NBS1 is a prostate cancer susceptibility gene. *Cancer Res* **64**, 1215–1219 (2004).
63. Friedman, J. M. Neurofibromatosis 1. in *GeneReviews®* (eds. Pagon, R. A. et al.) (University of Washington, Seattle, 1993).
64. Evans, D. G. R. *et al.* Cancer and Central Nervous System Tumor Surveillance in Pediatric Neurofibromatosis 1. *Clin Cancer Res* **23**, e46–e53 (2017).
65. Grolleman, J. E. *et al.* Mutational Signature Analysis Reveals NTHL1 Deficiency to Cause a Multi-tumor Phenotype. *Cancer Cell* **35**, 256–266.e5 (2019).
66. Yang, X. *et al.* Cancer Risks Associated With Germline PALB2 Pathogenic Variants: An International Study of 524 Families. *J Clin Oncol* **38**, 674–685 (2020).
67. Song, H. *et al.* Population-based targeted sequencing of 54 candidate genes identifies PALB2 as a susceptibility gene for high-grade serous ovarian cancer. *J Med Genet* **58**, 305–313 (2021).
68. Palles, C. *et al.* Germline mutations affecting the proofreading domains of POLE and POLD1 predispose to colorectal adenomas and carcinomas. *Nat. Genet.* **45**, 136–144 (2013).
69. Valle, L. *et al.* New insights into POLE and POLD1 germline mutations in familial colorectal cancer and polyposis. *Hum Mol Genet* **23**, 3506–3512 (2014).
70. Palles, C. *et al.* The clinical features of polymerase proof-reading associated polyposis (PPAP) and recommendations for patient management. *Fam Cancer* (2021) doi:10.1007/s10689-021-00256-y.
71. Logan, C. V. *et al.* DNA Polymerase Epsilon Deficiency Causes IMAGE Syndrome with Variable Immunodeficiency. *Am J Hum Genet* **103**, 1038–1044 (2018).
72. Henry, M.-L., Osborne, J. & Else, T. POT1 Tumor Predisposition. in *GeneReviews®* (eds. Adam, M. P. et al.) (University of Washington, Seattle, 1993).
73. Lim, T. L. *et al.* Germline POT1 variants can predispose to myeloid and lymphoid neoplasms. *Leukemia* **36**, 283–287 (2022).
74. Mester, J. & Charis, E. PTEN hamartoma tumor syndrome. *Handb Clin Neurol* **132**, 129–137 (2015).
75. Yang, X. *et al.* Ovarian and Breast Cancer Risks Associated With Pathogenic Variants in RAD51C and RAD51D. *J Natl Cancer Inst* **112**, 1242–1250 (2020).
76. Cybulski, C. *et al.* Germline RECQL mutations are associated with breast cancer susceptibility. *Nat. Genet.* **47**, 643–646 (2015).
77. Rohlin, A. *et al.* GREM1 and POLE variants in hereditary colorectal cancer syndromes. *Genes Chromosomes Cancer* **55**, 95–106 (2016).
78. Jaeger, E. *et al.* Hereditary mixed polyposis syndrome is caused by a 40-kb upstream duplication that leads to increased and ectopic expression of the BMP antagonist GREM1. *Nat. Genet.* **44**, 699–703 (2012).
79. Lieberman, S. *et al.* Features of Patients With Hereditary Mixed Polyposis Syndrome Caused by Duplication of GREM1 and Implications for Screening and Surveillance. *Gastroenterology* **152**, 1876–1880.e1 (2017).
80. Else, T., Greenberg, S. & Fishbein, L. Hereditary Paraganglioma-Pheochromocytoma Syndromes. in *GeneReviews®* (eds. Adam, M. P. et al.) (University of Washington, Seattle, 1993).
81. Grnøborg, S. *et al.* Leukoencephalopathy due to Complex II Deficiency and Bi-Allelic SDHB Mutations: Further Cases and Implications for Genetic Counselling. *JIMD Rep* **33**, 69–77 (2017).
82. Jackson, C. B. *et al.* Mutations in SDHD lead to autosomal recessive encephalomyopathy and isolated mitochondrial complex II deficiency. *J. Med. Genet.* **51**, 170–175 (2014).
83. Wain, K. E. *et al.* Appreciating the broad clinical features of SMAD4 mutation carriers: a multicenter chart review. *Genet. Med.* **16**, 588–593 (2014).
84. Jaspersion, K. & Burt, R. W. The Genetics of Colorectal Cancer. *Surg. Oncol. Clin. N. Am.* **24**, 683–703 (2015).
85. Tomlinson, I. P. & Houlston, R. S. Peutz-Jeghers syndrome. *J. Med. Genet.* **34**, 1007–1011 (1997).
86. Hearle, N. *et al.* Frequency and spectrum of cancers in the Peutz-Jeghers syndrome. *Clin. Cancer Res.* **12**, 3209–3215 (2006).
87. Giardiello, F. M. *et al.* Very high risk of cancer in familial Peutz-Jeghers syndrome. *Gastroenterology* **119**, 1447–1453 (2000).
88. Mai, P. L. *et al.* Risks of first and subsequent cancers among TP53 mutation carriers in the National Cancer Institute Li-Fraumeni syndrome cohort. *Cancer* **122**, 3673–3681 (2016).
89. Northrup, H., Krueger, D. A., & International Tuberous Sclerosis Complex Consensus Group. Tuberous sclerosis complex diagnostic criteria update: recommendations of the 2012 International Tuberous Sclerosis Complex Consensus Conference. *Pediatr. Neurol.* **49**, 243–254 (2013).
90. Krueger, D. A., Northrup, H., & International Tuberous Sclerosis Complex Consensus Group. Tuberous sclerosis complex surveillance and management: recommendations of the 2012 International Tuberous Sclerosis Complex Consensus Conference. *Pediatr. Neurol.* **49**, 255–265 (2013).
91. Au, K. S. *et al.* Genotype/phenotype correlation in 325 individuals referred for a diagnosis of tuberous sclerosis complex in the United States. *Genet. Med.* **9**, 88–100 (2007).
92. Salama, Y. *et al.* Comprehensive characterization of a Canadian cohort of von Hippel-Lindau disease patients. *Clin Genet* **96**, 461–467 (2019).
93. van Leeuwen, R. S., Ahmad, S., Links, T. P. & Giles, R. H. Von Hippel-Lindau Syndrome. in *GeneReviews®* (eds. Adam, M. P. et al.) (University of Washington, Seattle, 1993).